## Diffraction Grating

**BACKGROUND:** When coherent monochromatic light, such as that from a laser, passes through narrow slits an interference pattern is formed. A diffraction grating is composed of a large number of narrow evenly spaced slits. When laser light passes through the grating, a regular pattern of sharp bright maxima in the intensity of the light can be formed on a screen. The location of the mth maximum in the pattern is given by the relationship

$$m\lambda = d\sin(\theta_m), \qquad m = 1, 2, 3, \dots \tag{1}$$

where m is the order of the diffraction maximum, d is the separation between slits,  $\lambda$  is the wavelength of the light, and  $\theta_m$  is the angular displacement from the center of the zeroth order maximum (center of the pattern) to the center of the  $m^{\text{th}}$  order maximum.



Figure 9.1: Schematic arrangement of diffraction experiment.

If the linear distance between the center of the zeroth order maximum and the mth bright maximum is given by  $x_m$ , then the  $sin(\theta_m)$  is approximately given by

$$\sin(\theta_{\rm m}) = \frac{x_m}{L} \tag{2}$$

where L is the distance from the grating to the screen where the pattern is displayed and  $L \gg x_m$ . Combining Eqs. (1) and (2) yields

$$x_m = \frac{m\lambda L}{d}$$
, or  $\lambda = \frac{dx_m}{mL}$ , for  $m = 1, 2, 3...$  (3)

If the diffraction angle,  $\theta_m$ , is not small, then the wavelength must be calculated from the formula:

$$\lambda = \frac{r}{m} \sin\left[\tan^{-1}\left(\frac{x_m}{L}\right)\right] \,. \tag{4}$$

**OBJECT:** To study the characteristics of a diffraction grating and to measure the wavelength of the light from the HeNe laser.

**APPARATUS:** Diffraction Mosaic, meter stick, He Ne Laser ( $\lambda$ =632.8 nm).

**WARNING:** Never look directly into a laser beam.

Diffraction Grating:



Figure 9.2: The Diffraction Mosaic.

## PROCEDURE:

- 1. Arrange the laser and diffraction grating mosaic so that the mosaic is about 5m from a white wall in the laboratory. This wall will serve as the screen.
- 2. Place the laser behind the Mosaic so that the beam is incident normal to the 100 lines/mm grating (grating in upper left corner of the Mosaic).
- 3. Measure the separations between the zeroth order maximum and the 1st, 2nd and 3rd order maxima.
- 4. Record these data in the data table and determine the value for the wavelength,  $\lambda$ .
- 5. Compare your average calculated value of the wavelength to the given value for the HeNe laser.
- 6. Repeat steps 3 and 4 for the 300 lines/mm grating (top center of Mosaic).
- 7. Compare your average calculated value of the wavelength to the given value for the HeNe laser.
- 8. Repeat steps 3 and 4 for the 600 lines/mm grating (top right of Mosaic).
- 9. Compare your average calculated value of the wavelength to the given value for the HeNe laser.

Grating	d	L	Xm	tan(θm)=Xm/L	$\theta_m = \tan^{-1}(X_m/L)$	$\lambda$ =dsin $\theta_m$ /m
100 lines/mm	1mm/100		<b>X</b> 1		θ1 =	
			<b>X</b> 2		$\theta_2 =$	
			<b>X</b> 3		θ <sub>3</sub> =	
300 lines/mm	1mm/300		<b>X</b> 1		θ1 =	
			<b>X</b> 2		$\theta_2 =$	
			<b>X</b> 3		θ <sub>3</sub> =	
600 lines/mm	1mm/600		<b>X</b> 1		θ1 =	
			<b>X</b> 2		$\theta_2 =$	
			<b>X</b> 3		$\theta_3 =$	

## **Questions:**

- 1. Describe the differences observed in the diffraction patterns for the three gratings. 2. Show that Eq. (2) is (or is not) a valid approximation when  $L \gg x_m$ .